Europe / Middle-East / Africa

Microlife AG

Espenstrasse 139

9443 Widnau / Switzerland

Tel. +41 / 71 727 70 30

Fax +41 / 71 727 70 39

Email admin@microlife.ch

www.microlife.com

Asia

Microlife Corporation.

9F, 431, RuiGang Road, NeiHu

Taipei, 114, Taiwan, R.O.C.

Tel. 886 2 8797-1288

Fax 886 2 8797-1283

Email service@microlife.com.tw

www.microlife.com

North / Central / South America

Microlife USA, Inc.

424 Skinner Blvd., Suite C

Dunedin, FL 34698 / USA

Tel. +1 727 451 0484

Fax +1 727 451 0492

Email msa@microlifeusa.com

www.microlife.com

microlife BP 2BHO

Automatic Blood Pressure Monitor

Instruction Manual (1-11)

Tensiómetro automatico

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Table of contents

1. Introduction

- 1.1. Features
- 1.2. Important information about self-measurement

2. Important information on the subject of blood-pressure and its measurement

- 2.1. How does high/low blood-pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?

3. The various components of the blood-pressure monitor

4. Putting the blood-pressure monitor into operation

- 4.1. Inserting the batteries
- 4.2. Using a mains adapter
- 4.3. Tube connection

5. Carrying out a measurement

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Setting the inflation pressure
- 5.5. Measuring procedure
- 5.6. Discontinuing a measurement
- 5.7. Memory – displaying the last measurements

6. Error messages/malfunctions

7. Care and maintenance, recalibration

8. Guarantee

9. Technical specifications

10. www.microlife.com

1. Introduction

1.1. Features

The blood-pressure monitor is a fully automatic, digital blood-pressure measuring device for use on the upper arm, which enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring.

The device offers a very high and clinical tested measurement accuracy and has been designed to provide a maximum of user-friendliness.

Thanks to the large LCD display, both operational status and cuff pressure are visible for the entire duration of the measurement.

Before using, please read through this instruction manual carefully and then keep it in a safe place.

For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

Attention!

1.2. Important information about self-measurement

- Do not forget: **self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

Electromagnetic interference:

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

2. Important information on the subject of blood-pressure and its measurement

2.1. How does high/low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory centre, and adapted to the respective situation by way of feedback via the nervous system.

To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood-vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: During the «blood ejection» (Systole) the value is maximal (systolic blood-pressure value), at the end of the heart's «rest period» (Diastole) minimal (diastolic blood-pressure value).

The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2 Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood-pressure is over 140 mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

With blood-pressure values that are too low, i.e. systolic values under 105 mmHg and/or diastolic values under 60 mmHg, likewise, please consult your doctor.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.**

Table for classifying blood-pressure values (units mmHg) :

Range	Systolic Blood-pressure	Diastolic Blood-pressure	Measures
Hypotension	lower than 105	lower than 60	Consult your doctor
Normal range	between 105 and 140	between 60 and 90	Self-check
Limit value hypertension	between 140 and 160	between 90 and 95	Consult your doctor
Mild hypertension	higher than 160	between 95 and 105	Consult your doctor
Moderately serious hypertension	higher than 160	between 105 and 115	Consult your doctor
Serious hypertension	higher than 160	higher than 115	Consult your doctor urgently

Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». Please consult your doctor if you suspect that this might be the case.
- Correctly measured diastolic blood-pressure values above 120 mmHg require immediate medical treatment.

2.3. What can be done, if regular increased/low values are obtained?

- a) Please consult your doctor.
- b) Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result. Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.
- c) There are many different causes of the appearance of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take, not only for reducing a medically established high blood pressure, but also for prevention. These measures are part of your general way of life:

A) Eating habits

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

Follow consistently any medical instructions for treating previous illness such as:

- Diabetes (Diabetes mellitus)
- Fat metabolism disorder
- Gout

C) Habits

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)


D) Physical constitution

- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

3. The various components of the blood-pressure monitor

The illustration shows the blood-pressure monitor, consisting of:

a) Measuring unit


b) Upper arm cuff

Type M for arm circumference 22-32 cm

Type L for arm circumference 32-42 cm


(available as special accessory)


4. Putting the blood-pressure monitor into operation

4.1. Inserting the batteries

After you have unpacked your device, first insert the batteries. The battery compartment is located on the underside of the device (see illustration).

- Remove cover as illustrated
- Insert the batteries (4 x size AA 1,5 V), thereby observing the indicated polarity.
- If the battery warning appears in the display, the batteries are empty and must be replaced by new ones.

Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AA» Long-Life or Alkaline 1.5V Batteries. The use of 1.2V Accumulators is not recommended.
- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.

Functional check:


Hold on the O/I button down to test all the display elements. When functioning correctly all segments must appear.

4.2. Using a mains adapter (special accessory)

It is possible to operate this blood-pressure instrument through a mains adapter (output 6V DC / 600 mA with DIN plug). Make certain that you use a mains adapter which fulfils the legal requirements (CE sign on the type plate).

- Plug the DIN plug into the socket at the back of the instrument.
- Plug the mains adapter into a 230 V power socket.
- Test that power is available by pressing the O/I button.


Note:

- No power is taken from the batteries while the mains adapter is connected to the instrument.
- If the mains voltage is interrupted during the measurement (e.g. by removal of the mains adapter from the mains socket), the instrument must be reset by removing the DIN plug from its socket.
- Please consult a specialist dealer if you have questions relating to the mains adapter.

4.3. Tube connection

Insert the cuff tube into the opening provided on the left side of the instrument, as shown in the diagram.

5. Carrying out a measurement

5.1. Before the measurement

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about 5 minutes before the measurement.
- Remove any garment that fits closely to your upper arm.
- Measure always on the same arm (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.


5.2. Common sources of error

Note: Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood-pressure will be measured! (Each 15 cm difference in height results in a measurement error of 10mmHg!)
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is of extraordinary importance. The cuff size is dependent upon the circumference of the arm (measured in the centre). The permissible range is printed on the cuff. If this is not suitable for your use, please contact your dealer. **Note:** Only use clinically approved Original-Cuffs!
- A loose cuff or a sideways protruding air-pocket causes false measurement values.

5.3. Fitting the cuff

- a) Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.
- b) Make certain that the lower edge of the cuff lies approx. 2 to 3 cm above the elbow and that the rubber tube leaves the cuff on the inner side of the arm.


- Important!** The mark (ca. 3 cm long bar) must lie over the artery which runs down the inner side of the arm.
- c) Tighten the free end of the cuff and close the cuff with the closer.
 - d) There must be no free space between the arm and the cuff as this would influence the result. Clothing must not restrict the arm. Any piece of clothing which does (e.g. a pullover) must be taken off.

- e) Secure the cuff with the Velcro closer in such a way that it lies comfortably and is not too tight. Lay the arm on the table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked.

- f) Remain seated quietly for 2 minutes before you begin the measurement.

Comment:




If it is not possible to fit the cuff to the left arm, it can also be placed on the right one. However all measurements should be made using the same arm

5.4. Setting the inflation pressure

Program the inflation pressure to be approximately 40 mmHg higher than your estimated systolic blood pressure. You can choose between 8 initial values: 140, 160, 180, 200, 220, 240, 260 and 280 mmHg. If you do not know your systolic blood pressure, select a value between 160 and 200. The instrument is pre-set to 140 mmHg. If this initial value is sufficient for you, omit the following step and start immediately with the measurement.

Set the initial pressure as follows:

- Switch on the instrument by pressing the O/I button. This picture appears on the LCD display. (It remains there as long as you hold the button pressed down).
- After the button is released, the display flashes as shown and the instrument beeps. This indicates that the instrument is automatically adjusting itself to the pressure of the surroundings.
- After this equalization is completed, zero is displayed. Finally a beep is emitted to indicate that the instrument is ready to take measurements.
- On pressing the SET button the instrument displays the pre-set initial pressure for approx. 3 seconds.
- In order to set the initial pressure individually, press the SET button briefly for each level until the required initial pressure is reached. Each button press raises the initial pressure by 20 mmHg. After reaching 280 mmHg the display starts again at 120 mmHg.
- After the SET button has been released, the instrument is set to the displayed value and the setting up procedure is complete. As long as the batteries are not removed, the value is stored in the instrument even after it has been switched off.


Comment:

If the initial pressure selected was too low it will be corrected automatically during the measurement by additional inflation. This results in a slightly longer measuring procedure which the patient may find uncomfortable. The result of the measurement is normally not affected.

5.5. Measuring procedure

After the cuff has been appropriately positioned and the setting of the pressure pre-selection has been undertaken, the measurement can begin:


- Press O/I button to put the device into operation.
- Press the START button. The pump begins to inflate the cuff.
The rising pressure in the cuff is shown in the display.
- After the inflation pressure has been reached, the pump stops and the pressure gradually falls.
The cuff pressure are displayed.
- When the instrument detects a pulse, the heart symbol in the display starts to flash and a beep is heard for every heart beat.
- A longer beep is sounded when the measurement has been completed. The systolic, diastolic blood pressures and the pulse rate of the patient now appear in the display.
- The results of the measurement remain on the display until you switch off the instrument. If no button is pressed for a period of 5 minutes the instrument switches itself off in order to preserve the batteries.

5.6. Discontinuing a measurement


If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the O/I button can simply be pressed at any time.

The instrument then immediately reduces the cuff pressure.

5.7. Memory – displaying the last measurement

The measured results are stored in the instrument until a new measurement is carried out or the batteries are removed. You can call up the stored values at any time by pressing the MEMORY button.

6. Error messages/malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.
(example: error no. 3)

Error No.	Possible cause(s)
ERR 1	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. The tube may have become unplugged after the systolic blood pressure was measured. Further possible cause: No pulse could be detected
ERR 2	Unnatural pressure impulses impair the measurement result. Possible cause: The arm was moved during the measurement (artefact).
ERR 3	Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.
HI	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
LO	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

* Please consult your doctor, if this or any other problem occurs repeatedly.

Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched on although the batteries are in place.	<ol style="list-style-type: none"> 1. Check whether the batteries are installed with the correct polarity and correct if necessary. 2. If the display is unusual, remove the batteries and then exchange them for new.
The pressure does not rise although the pump is running.	<ul style="list-style-type: none"> • Check the connection of the cuff tube and connect properly if necessary.
The instrument frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none"> 1. Fit the cuff correctly on the arm. 2. Before starting the measurement make sure that the cuff is not fitted too tightly and that there is no rolled-up sleeve exerting pressure on the arm above the measuring position. Take off articles of clothing if necessary. 3. Measure the blood pressure again in complete peace and quiet.
Every measurement produces a different value although the instrument functions normally and the values displayed are normal.	<ul style="list-style-type: none"> • Please read the following information and the points listed under «Common sources of error». Repeat the measurement.
The blood pressure values measured differ from those measured by the doctor.	<ul style="list-style-type: none"> • Record the daily development of the values and consult your doctor about them.

Further information

The level of blood-pressure is subject to fluctuations even with healthy people. Important thereby is, that **comparable measurements always require the same conditions (Quiet conditions)!**

You must consult your specialist dealer or chemist if there are technical problems with the blood-pressure instrument. **Never attempt to repair the instrument yourself!**

Any unauthorised opening of the instrument invalidates all guarantee claims!

7. Care and maintenance, recalibration

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. **The cuff must not be washed!**
- d) Handle the tube carefully. Do not pull on it. Likewise, this is not to be buckled or laid over sharp edges.
- e) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the device!** Otherwise the manufacturer calibration becomes invalid!

Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display **every 2 years.**

Your specialist dealer would be pleased to provide more extensive information about this.

8. Guarantee

The blood-pressure monitor is **guaranteed for 3 years** from date of purchase. This guarantee includes the instrument and the cuff. The guarantee does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument by third parties.

The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

Name and company address of the responsible dealer:

9. Technical specifications

Weight:	575 g (with batteries + cuff size M)
Size:	106 (W) x 162 (L) x 61.5 (H) mm
Storage temperature:	-5 to +50 °C
Humidity:	15 to 90% relative humidity maximum
Operation temperature:	10 to 40 °C
Display:	LCD-Display (L iquid C rystal D isplay)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0–299 mmHg
Memory:	Storing the last measurements automatically
Measuring resolution:	1 mmHg
Accuracy:	
Pressure:	within ± 3 mmHg
Pulse:	± 5 % of the reading
Power source:	a) 4 x 1.5 V Batteries; size AA b) Mains adapter DC 6V, 600 mA (optional)
Cuffs:	type M for arm circumference 22–32 cm or type L for arm circumference 32–42 cm (optional)
Reference to standards:	EU Directives 93/42/EEC; NIBP requirements: EN 1060-1 /-3 /-4, ANSI / AAMI SP10

Technical alterations reserved!

10. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com.

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