*Read the instructions carefully before using this device.*

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***microlife***<sup>®</sup>

# Wrist Watch Blood Pressure Monitor

## Instruction Manual

### Important Safety Instructions

Type BF applied part

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.

---

**1. Introduction**

- 1.1. Features
- 1.2. Important information about self-measurement

---

**2. Important information on the subject of blood-pressure and its measurement**

- 2.1. How does high/low blood-pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?

---

**3. The various components of the blood-pressure monitor**

---

**4. Putting the blood-pressure monitor into operation**

- 4.1. Inserting the batteries
- 4.2. Reading the set date
- 4.3. Setting the time and date

---

**5. Carrying out a measurement**

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Discontinuing a measurement
- 5.6. Memory – storage and recall of the measurements

---

**6. Error messages/malfunctions**

---

**7. Care and maintenance, recalibration**

---

**8. Guarantee**

---

**9. Technical specifications**

---

**10. [www.microlife.com](http://www.microlife.com)**

---

## 1. Introduction

---

### 1.1. Features

The blood-pressure monitor (with integrated time/date display) is a fully automatic, digital blood-pressure measuring device for use on the wrist, which enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measuring.

The device offers a very high and clinical tested measurement accuracy and has been designed to provide a maximum of user-friendliness.

Before using, please read through this instruction manual carefully and then keep it in a safe place. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

### Attention!

#### 1.2. Important information about self-measurement

- Do not forget: self-measurement means **control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.
- The pulse display is not suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

#### Electromagnetic interference:

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

## 2. Important information on the subject of blood-pressure and its measurement

---

### 2.1. How does high/low blood-pressure arise?

The level of blood-pressure is determined in a part of the brain, the so-called circulatory centre, and adapted to the respective situation by way of feedback via the nervous system. To adjust the blood-pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. The latter is effected by way of fine muscles in the blood-vessel walls.

The level of arterial blood-pressure changes periodically during the heart activity: During the «blood ejection» (Systole) the value is maximal (systolic blood-pressure value), at the end of the

heart's «rest period» (Diastole) minimal (diastolic blood-pressure value). The blood-pressure values must lie within certain normal ranges in order to prevent particular diseases.

## 2.2 Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood-pressure is over 140 mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

With blood-pressure values that are too low, i.e. systolic values under 100 mmHg and/or diastolic values under 60 mmHg, likewise, please consult your doctor.

Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to alter independently the drug doses prescribed by your doctor.**

Table for classifying blood-pressure values (units mmHg) according to World Health Organization:

Range	Systolic Blood-pressure	Diastolic Blood-pressure	Measures
Hypotension	lower than 100	lower than 60	Consult your doctor
Normal range	between 100 and 140	between 60 and 90	Self-check
Mild hypertension	between 140 and 160	between 90 and 100	Consult your doctor
Moderate hypertension	between 160 and 180	between 100 and 110	Consult your doctor
Severe hypertension	higher than 180	higher than 110	Consult your doctor immediately

### Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». Please consult your doctor if you suspect that this might be the case.
- Correctly measured diastolic blood-pressure values above 120 mmHg require **immediate medical treatment.**

### **2.3. What can be done, if regular increased/low values are obtained?**

- a) Please consult your doctor.
- b) Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result. Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.
- c) There are many different causes of the appearance of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take, not only for reducing a medically established high blood pressure, but also for prevention. These measures are part of your general way of life:

#### **A) Eating habits**

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

#### **B) Previous illnesses**

Follow consistently any medical instructions for treating previous illness such as:

- Diabetes (Diabetes mellitus)
- Fat metabolism disorder
- Gout

#### **C) Habits**

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)


#### **D) Physical constitution:**

- After a preliminary medical examination, **do regular sport**.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

### 3. The various components of the blood-pressure monitor

---

The illustration shows the blood-pressure monitor, consisting of:


Wrist cuff Type WC1 13.5–19.5 cm, for wrist circumference of 13.5 to 19.5 cm

### 4. Putting the blood-pressure monitor into operation

---

#### 4.1. Inserting the batteries

After you have unpacked your device, first insert the batteries. The battery compartment is located on the left lateral side of the device (see illustration).

- Remove cover as illustrated
- Insert the batteries (2 x size AAA 1.5V), thereby observing the indicated polarity.
- If the battery warning appears in the display, the batteries are empty and must be replaced by new ones.

#### Attention!


- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AAA» Long-Life or Alkaline 1.5V Batteries. The use of 1.2V Accumulators is not recommended.

- If the blood-pressure monitor is left unused for long periods, please remove the batteries from the device.

**Functional check:** Hold on the O/I button down to test all the display elements. When functioning correctly all segments must appear.

#### 4.2. Reading the set date

Please press the Time button.

#### 4.3. Setting the time and date

- 1) After the new batteries are fitted, the year number flashes in the display. You can set the year by pressing the M-button. To confirm and then set the month, press the Time button.
- 2) You can now set the month using the M-button. Press the time button to confirm and then set the day.
- 3) Please follow the instructions above to set the day, hour and minutes.
- 4) Once you have set the minutes and pressed the Time button, the date and time are set and the time is displayed.
- 5) If you want to change the date and time, press and hold the Time button down for approx. 3 seconds until the year number starts to flash. Now you can enter the new values as described above.

### 5. Carrying out a measurement

---

#### 5.1. Before the measurement

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about 5 minutes before the measurement.
- Measure always on the same wrist (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

#### 5.2. Common sources of error


**Note:** Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.

- If the wrist artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood-pressure will be measured! (Each 15 cm difference in height results in a measurement error of 10 mmHg!)
- A loose cuff causes false measurement values.

### 5.3. Fitting the cuff

- Remove all eventual objects and jewellery (e.g. wrist-watch) from the wrist in question. Draw the cuff over the wrist.
- The distance between the cuff and the hand should be approx. 10 mm.
- Secure the cuff with the fastener, so that it lies comfortably and **not too tight**, whereby no space should remain between the cuff and the wrist.
- Lay the arm on a table, with the palm upwards. Support the arm a little with a rest (cushion), so that the cuff rests at about the same height as the heart. Take care, that the cuff lies free.

### 5.4. Measuring procedure

After the cuff has been appropriately positioned and the setting of the pressure pre-selection has been undertaken, the measurement can begin:

- Press O/I button to put the device into operation
- Press the START button – the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.
- After reaching the inflation pressure, the pump stops and the pressure slowly falls away. The cuff-pressure (large characters) is displayed during the measurement. When the device has detected the pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.

START

- d) When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values as well as the pulse frequency are now displayed.

Example (Fig.): Systole 118, Diastole 73, Pulse 75


The measurement results are displayed, until you switch the device off. If no button is pressed for 5 minutes, the device switches automatically off, to save the batteries.

### 5.5. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the O/I button can be pressed at any time.

The device then immediately lowers the cuff-pressure automatically.

### 5.6. Memory – storage and recall of the measurements


The blood pressure monitor automatically stores each of the last 30 measurement values. By pressing the M-button, the last measurement (MR1) as well as the further last 29 measurements (MR2, MR3, ..., MR30) can be displayed one after the other.

(MR1: Values of the last measurement) (MR2-MR30: Values of the measurement before MR1)

## 6. Error messages/malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).

Error No.	Possible cause(s)
ERR 1	No pulse has been detected.
ERR 2	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artefact).
ERR 3	The inflation of the cuff takes too long. The cuff is not correctly seated.
ERR 5	The measured reading indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully.*
HI	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
LO	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

\* Please consult your doctor, if this or any other problem occurs repeatedly.

**Other possible malfunctions and their elimination** – If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched on although the batteries are in place.	<ol style="list-style-type: none"> <li>1. Check whether the batteries are installed with the correct polarity and correct if necessary.</li> <li>2. If the display is unusual, re-insert batteries or exchange them.</li> </ol>
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none"> <li>1. Check the positioning of the cuff.</li> <li>2. Measure the blood pressure again in complete peace and quiet under observance of the details made under point 5.</li> </ol>
Every measurement results in different values although the device functions normally and the values displayed are normal.	<ul style="list-style-type: none"> <li>• Please read the following information and the points listed under «Common sources of error». Repeat the measurement.</li> </ul>
The blood pressure values measured differ from those measured by the doctor.	<ul style="list-style-type: none"> <li>• Record the daily development of the values and consult your doctor about them.</li> </ul>

### Further Information

The level of blood pressure is subject to fluctuations even with healthy people. Important thereby is, that comparable **measurements always require the same conditions (rest conditions)!**

You must consult your specialist dealer or chemist if there are technical problems with the blood pressure instrument. **Never attempt to repair the instrument yourself!**

Any unauthorised opening of the instrument invalidates all guarantee claims!

## 7. Care and maintenance, recalibration

---

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!
- d) **Do not drop** the instrument or treat it roughly in any way. Avoid strong vibrations.
- e) Never open the device! Otherwise the manufacturer calibration becomes invalid!

### Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display **every 2 years**. Your specialist dealer would be pleased to provide more extensive information about this.

## 8. Guarantee

---

The blood-pressure monitor **guaranteed for 3 years** from date of purchase. This guarantee includes the instrument and the cuff. The guarantee does not apply to damage caused by improper handling, accidents, not following the operating instructions or alterations made to the instrument by third parties.

The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

Name and company address of the responsible dealer:

## 9. Technical specifications

---

<b>Weight:</b>	145g (with batteries + cuff)
<b>Size:</b>	80 x 70 x 75 mm (including cuff)
<b>Storage temperature:</b>	-20 - +55 °C; -4 - +131 °F; 15 - 95 % relative maximum humidity
<b>Operation temperature:</b>	10 - 40 °C; 50 - 104 °F; 15 - 95 % relative maximum humidity
<b>Display:</b>	LCD-Display (Liquid Crystal Display)
<b>Measuring method:</b>	oscillometric
<b>Pressure sensor:</b>	capacitive
<b>Measuring range:</b>	
<b>SYS/DIA:</b>	30 to 280 mmHg
<b>Pulse:</b>	40 to 200 per minute
<b>Cuff pressure display range:</b>	0-299 mmHg
<b>Memory:</b>	Storing the last 30 measurements automatically
<b>Measuring resolution:</b>	1 mmHg
<b>Pressure accuracy:</b>	Pressure within $\pm 3$ mmHg Pulse $\pm 5$ % of the reading
<b>Power source:</b>	2 x 1.5 V Batteries; size AAA
<b>Accessories:</b>	storage case
<b>Reference to standards:</b>	EN 1060-1 /-3 /-4; IEC 60601-1; IEC 60601-1-2 (EMC)

This device complies with the requirements of the Medical Device Directive 93/42/EEC.  
Technical alterations reserved.

## 10. [www.microlife.com](http://www.microlife.com)

---

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at [www.microlife.com](http://www.microlife.com).

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