Bedienungsanleitung/Garantie

Gebruiksaanwijzing • Mode d'emploi • Instrucciones de servicio
Manual de instruções • Istruzioni per l'uso • Bruksanvisning
Instruction Manual • Instrukcja obsługi/Gwarancja
Návod k použití • Használati utasítás • Руководство по эксплуатации

Munddusche

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General Safety Instructions

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing. If you give this device to other people, please also pass on the operating instructions.

- The appliance is designed exclusively for private use and for the envisaged purpose. This appliance is not fit for commercial use. Do not use it outdoors. Keep it away from sources of heat, direct sunlight, humidity (never dip it into any liquid) and sharp edges. Do not use the appliance with wet hands. If the appliance is humid or wet, unplug it immediately.
- When cleaning or putting it away, switch off the appliance and always pull out the plug from the socket (pull the plug itself, not the lead) if the appliance is not being used and remove the attached accessories.
- Do **not** operate the machine without supervision. If you leave the room you should always turn the device off. Remove the plug from the socket.
- Check the appliance and the cable for damage on a regular basis. Do not use the appliance if it is damaged.
- Do not try to repair the appliance on your own. Always contact an authorized technician. To avoid the exposure to danger, always have a faulty cable be replaced only by the manufacturer, by our customer service or by a qualified person and with a cable of the same type.
- Use only original spare parts.
- Pay careful attention to the following "Special Safety Instructions".

Children and Frail Individuals

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.

WARNING!

Caution! Do not allow small children to play with the foil as there is a **danger of suffocation!**

- In order to protect children or frail individuals from the hazards of electrical equipment, please ensure that this device is only used under supervision. This device is not a toy. Do **not** allow small children to play with it.

Symbols in these Instructions for Use

Important information for your safety is specially marked. It is essential to comply with these instructions in order to avoid accidents and prevent damage to the machine:

WARNING:

This warns you of dangers to your health and indicates possible injury risks.

CAUTION:

This refers to possible hazards to the machine or other objects.

NOTE: This highlights tips and information.

Special Safety Instructions

- Please use the device only for cleaning the mouth.
- Please empty the water container after use.
- The device should not be immersed in water.
- Ensure that the water hose (the connection between the device and the water jet) is not squeezed
- If you wear a brace, please ask your dentist if it is possible to use a dental water jet. We do not accept liability for any damage caused.
- If you have any form of disease in the mouth or throat, please ask your doctor whether the use of a dental water jet is advisable.

Overview of the Components

- 1 Adjustment ring (wet jet or soft jet)
- 2 Single jet nozzle
- 3 Eject button for single jet nozzle
- 4 Water stop
- 5 Water hose
- 6 On/off switch
- 7 Water pressure regulator
- 8 Water container

Wall Installation

WARNING OF ELECTRIC SHOCK!

- The base of the device is not protected against immersion in water.
- For your safety the charging station should therefore be firmly installed.

Install the station on the wall. Please remember that there must be an electric socket close by.

- Please measure the distance between the indentations on the rear and drill two holes for the rawlplugs the same distance apart.

CAUTION:

Please ensure beforehand that you will not be damaging any pipes or electric wires concealed in the wall!

Unpacking the Device

- Take the device out of its packaging.
- Remove all packaging materials, such as foils, filler material, cable brackets and cardboard packaging.
- Check the scope of delivery.

NOTE:

There may still be some dust or production residues on the device. We recommend that you clean the device as described under „Cleaning“.

Before Initial Use

- Fill the water container with water.
- Turn the on/off switch (6) on.
- Turn the regulator to the maximum setting. Then switch on the device on the handle.
- With the adjustment ring (1) you can set the type of water jet. You have the choice between a soft jet and a massage jet.
- Allow the water to flow directly into a wash basin in order to remove any residues.
- Rinse the nozzle attachments briefly under running water.

Nozzles

The scope of delivery contains 4 nozzles with different coloured rings. These can be allocated to different users.

Operation

- Before using the water jet you should clean your teeth with a toothbrush.
- Attach a nozzle on the handle in such a way that the aperture of the nozzle is pointing at the water stop button. Ensure that the nozzle is fixed firmly in position.
- Fill the water container with water which is hand warm. Please do not fill in any hot water.
- In order to clean the nozzles, allow water to flow through the nozzle head (at maximum water pressure) into the wash basin before and after each time the device is used.
- Then use the water pressure regulator to select the desired water pressure. Be-gin on the smallest setting, especially if you have problems with bleeding of the gums or when the water jet is used by children.
- Make sure that you hold your head above the wash basin. Pass the nozzle along the gums and teeth.
- After use, please turn the water pressure regulator back to position „Min“.
- Turn the device off on the handle.
- Please clean the device as described under „Cleaning and Care“.

Cleaning and Care

- Before cleaning the water jet, remove the mains plug from the socket. The housing and mains supply must not be cleaned or immersed in water. The housing should only be wiped off with a moist soft cloth.
- Pour away the remaining water from the water container.

- After use, always remove the nozzles from the handle and clean them thoroughly under warm running water, before drying them off and returning them to the storage container.
- Drive the handle before inserting it into the water jet.
- When water additives are used, rinse out the water container with clear water. Then fill the water container with 100 ml of clear water and pass the water through the device into the wash basin.

Technical Data

Model:.....MD 5503
 Power supply:.....220-240 V, 50 Hz
 Power consumption:.....15 W
 Protection class:.....II
 Net weight:.....0,98 kg

This device has been tested according to all relevant current CE guidelines, such as electromagnetic compatibility and low voltage directives, and has been constructed in accordance with the latest safety regulations.

Subject to technical changes without prior notice!

Meaning of the “Dustbin” Symbol

Protect our environment: do not dispose of electrical equipment in the domestic waste.

Please return any electrical equipment that you will no longer use to the collection points provided for their disposal.

This helps avoid the potential effects of incorrect disposal on the environment and human health.

This will contribute to the recycling and other forms of reutilisation of electrical and electronic equipment.

Information concerning where the equipment can be disposed of can be obtained from your local authority.

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