**Oral-B®**

**Vitality™**

D 12513

D 12013

powered by **BRAUN**

Type 3709

|            |                                |
|------------|--------------------------------|
| <b>UK</b>  | 0 800 73 11 792                |
| <b>IRL</b> | 1 800 509 448                  |
| <b>F</b>   | 0 810 365 855                  |
| <b>B</b>   | 0 800 14 592                   |
| <b>LU</b>  | 4 00 50 51                     |
| <b>PL</b>  | 0 801 127 286<br>0 801 1 BRAUN |
| <b>CZ</b>  | 221 804 335                    |
| <b>SK</b>  | 02/5710 1135                   |
| <b>H</b>   | 06-1-451-1256                  |
| <b>TR</b>  | 0 212 473 75 85                |
| <b>RUS</b> | 8-800-200-20-20                |
| <b>UA</b>  | +38 044 428 65 05              |

**Internet:**









[www.oralb.com](http://www.oralb.com)

[www.braun.com](http://www.braun.com)

[www.service.braun.com](http://www.service.braun.com)

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3-713-100/00/I-08/M  
GB/F/PL/CZ/SK/H/TR/RUS/UA/Arab

## English

Please read the use instructions carefully before first use.

### Important

- Periodically check the cord for damage. A damaged or non-functioning unit shall no longer be used. Take the charging unit (D) to an Oral-B Braun Service Centre.
- This product is not intended for use by children under age three.
- This appliance is not intended for use by children or persons with reduced physical or mental capabilities, unless they are given supervision by a person responsible for their safety. In general, we recommend that you keep the appliance out of reach of children.
- If the product is dropped, the brush head should be replaced before the next use even if no damage is visible.

## Description

- A Brush head (depending on model) – interchangeable
- B On/off switch
- C Handle
- D Charging unit

## Charging

- Plug the charging unit (D) into an electrical outlet. Place the handle on the charging unit. Handle must be switched off.
- A full charge takes 16 hours and will provide an operation time of approximately 20 minutes.
- For everyday use, the handle can be stored on the plugged-in charging unit to maintain the toothbrush at full power. Overcharging is impossible.

### Battery maintenance

To maintain the maximum capacity of the rechargeable battery, you must discharge it at least once every 6 months. To do this, unplug the charging unit and allow the battery to run down with normal use.

## Using the toothbrush

- Wet brush head and apply toothpaste.
- Place brush head on teeth, then switch on.
- Slowly guide the brush head from tooth to tooth, first the outsides, then the insides, finally the chewing surfaces. Do not press too hard or scrub.
- After two minutes (brushing time recommended by dentists) the timer (D 12513 model only) will signal by a short stutter.
- During the first days of use, gums may bleed slightly. Should bleeding persist after 2 weeks, contact your dentist.
- Replace brush head every 3 months, when bristles are worn (1c) or when colour on Indicator bristles (1b) has faded half-way down the bristles.

## Cleaning

After use, rinse the brush head under running water (2). Remove the brush head and rinse brush head and handle separately (3), then dry.

Subject to change without notice.

This product conforms to the European Directives EMC 2004/108/EC and Low Voltage 2006/95/EC.

This product contains rechargeable batteries. In the interest of protecting the environment, please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at an Oral-B Braun Service Centre or at appropriate collection points provided in your country.

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# English

## Guarantee

We grant 2 years limited guarantee on the product commencing on the date of purchase. Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose. This guarantee extends to every country where this appliance is supplied by Braun or its appointed distributor.

This guarantee does not cover: damage due to improper use, normal wear or use, especially regarding the brush heads, as well as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorised persons and if original Braun parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Oral-B Braun Customer Service Centre.

## For UK only:

This guarantee in no way affects your statutory rights.

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